

TITANIUM COLORS



PAPIER KANN DIE BRILLANZ DER FARBEN NUR ANNÄHERND DARSTELLEN, GERINGE ABWEICHUNGEN SIND MÖGLICH.
PAPER CAN ONLY ROUGHLY CAPTURE THE BRILLIANCE OF OUR COLORS, THEREFORE SMALL DEVIATIONS ARE POSSIBLE.

VERGLASUNGSHINWEIS

Zur Glasrandung bitte die Stützscheiben abtasten.

Größe des Korrektionsglases im Vergleich zur Stützscheibe: + 0.1 mm

M19: Konektor zur Glasrandung entfernen.

1 Nutbreite: 0.8 mm.

2 Nuttiefe: M11 & M19: 0.45 mm, M14 & M20: 0.4 mm

3 Mindestranddicke Korrektionsglas: 2.0 mm

Bitte beachten Sie, dass es sich bei allen angegebenen Maßen um Endmaße handelt!

Bitte messen Sie nach!

Bitte beachten Sie unsere Montageanleitung auf www.meyer-eyewear.de/downloads

GLAZING NOTICE

Please trace the demo lens to start edging process.

Size of prescription lens compared to demo lens: +0.1 mm

M19: please remove the Konektor for glazing.

1 Groove width: 0.8 mm.

2 Groove depth: M11 & M19: 0.45 mm, M14 & M20: 0.4 mm

3 Minimum edge thickness: 2.0 mm

Please note that all given dimensions are final dimensions!

Please measure!

Please pay attention to our glazing instruction at www.meyer-eyewear.com/downloads

MEYER BRILLENMANUFAKTUR GMBH

URSULINENSTRASSE 63A 66111

SAARBRÜCKEN - GERMANY

T +49 (0) 681 96 86 641

F +49 (0) 681 96 86 647

